

l'acn SCIENTIFIC LABORATORIES

Via 25 Aprile 9/13 - 20023 CERRO MAGGIORE (MI) ITALY

Tel. +39 0331 420303

Fax. +39 0331 420153

Web: www.acn.it E-mail: acn@acn.it

Quality Control Kit for radiopharmaceutical preparations

l'acn chromatography kits allow to test the radiochemical purity of specific Tc^{99m} -labeled radiopharmaceuticals. The strips dimensions can be sized according to the instructions contained in the liftlet of the radiopharmaceutical or in accordance with the customer's needs (ITLC-SG, TLC-SG, ITLC-SA, Whatman, TLC Aluminum Oxide, ...). Some solvents must be purchased separately due to hazardous material shipping restrictions.

The kits include:

- No. 1 (or 2) cylinder for the development
- No. 1 forcep
- Disposable pipette to collect the reagents
- Supports for chromatography
 - a thin layer (ITLC-SG/SA or TLC-SG)
 - chromatography paper
 - chromatographic column
 - membrane filter



The **Aluminum Breakthru Kit** provides a rapid, easy and inexpensive way to test aqueous solutions, particularly pertechnetate generator eluate, for trace quantities of aluminum